

Redox Titration Lab

Name _____ Partner _____ Date _____ Score _____

Purpose:

- To determine the concentration of an unknown $\text{KMnO}_{4(\text{aq})}$ solution.
- To develop skill with a burette:
- To understand the properties of a primary standard
- To experience colour change as a visible marker for the endpoint of a reaction

Materials: funnel, volumetric flask, beaker, Erlenmeyer flask, burette, burette brush, pipette, $\text{H}_2\text{SO}_{4(\text{aq})}$, $\text{FeSO}_4 \cdot 7\text{H}_2\text{O}_{(\text{aq})}$, $\text{KMnO}_{4(\text{aq})}$

- Prelab: The $\text{KMnO}_{4(\text{aq})}$ will be the titrant. Write a balanced redox reaction. Be sure to include voltage

-
- ❖ Color of Excess reagent _____
 - ❖ Color of limiting reagent _____
 - ❖ Color at equivalence point _____
 - ❖ Color at end point _____

- Calculate the mass of iron (II) sulfate hepta hydrate that is required to make 0.10 L of 0.10 mol/L solution

Procedure:

- ❖ Make the solution of iron (II) sulfate.
 - o Weigh out the correct mass of iron (II) sulfate.
 - o Dissolve it in 50 mL of the 5.0 mol/L acid solution provided. Use a BEAKER for this step
 - o Transfer the solution to a volumetric flask. Do sufficient rinsing with DISTILLED water
 - o Fill the volumetric flask to the 100 mL mark. Use distilled water and an eye dropper.
 - o Stopper and invert several times.
 - o Pour the solution out into a clean, dry beaker.
- ❖ Pipette 10 mL of the $\text{FeSO}_{4(\text{aq})}$ from the beaker into an Erlenmeyer flask. Be sure to place a white piece of paper under the flask. This will maximize the colors in the flask.
- ❖ Clean and prepare the burette for the excess reagent ($\text{KMnO}_{4(\text{aq})}$) that is provided. Be sure to use a funnel. Record the initial volume of the burette.
- ❖ Titrate the iron (II) sulfate solution with the KMnO_4 solution until the end point is reached. Record the final volume of the burette.
- ❖ Keep the first trial for color comparison.
- ❖ Do successive trials of this lab until you get three consistent volumes of titrant. (must be within 0.1 mL) Keep the trials and compare the colors at the end point.

Observations:

Make a clear and well labeled table to show your observations.

Analysis of DATA

- Calculate concentration of the $\text{KMnO}_{4(\text{aq})}$

- Why is acidified $\text{KMnO}_{4(\text{aq})}$ a poor primary standard solution?